## NOTICE/IFU :

### Trousse chirurgicale (Surgical set) IDCam-IDBio, IDAII, RBS3.

## F Trousse EN Set

REF

RBS3, TROUSSEIDALL, TROUSSEIDCAM

## F

### Recommandations de forage

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### Drilling recommendations

#### Maxilla :

- 1) 1st drilling Ø2 mm with irrigation 650 rpm.
- 2) 2nd and 3rd drilling without irrigation at 150 rpm. With irrigation at 650 rpm (with the RBS® drill) or at 1500-2000 rpm (with the 4 blades drill).
- 3) Rinse abundantly the socket with tobramycin, diluted with physiological serum before implant setting.

#### Mandible :

- 1) 1st drilling Ø2 mm with irrigation 650 rpm.
- 2) 2nd and 3rd drilling and following with irrigation at 650 rpm.
- 3) Rinse abundantly the socket with tobramycin, diluted with physiological serum before implant setting.

### Implant setting

- 1) Place the implant at 1.5 mm subcrestally whatever the implant line.
- 2) Screw the closing and healing caps at 5N.cm by hand. **Expert advice :** soak the implant, the closing cap and the healing cap in the tobramycin antibiotic (75 mg) diluted in 20 cl of physiological serum before placing it.

### Sterilization and cleaning of surgical drill sets (manual cleaning)

- 1) Use **powder-free** gloves, also for the assistant.
  - 2) Soak the surgical set with the instruments in a decontamination tank during 15 min.
  - 3) Remove all the instruments from the set.
  - 4) Clean each instrument with a brass brush (brass only) in order to remove any bone fragment.
  - 5) Put the instruments + the empty base plate in an ultrasonic bath containing a **powdered decontamination product** for 15 minutes.
  - 6) Change the gloves (power-free).
  - 7) Rinse both instruments and base plate in demineralized water **during 5 minutes. Renew water five times.**
  - 8) Dry all parts.
  - 9) Put the instruments back on the base plate while wearing **powder-free gloves.**
  - 10) Sterilize by autoclaving at 134°C 18 minutes according to the instructions of European Pharmacopoeia 8th edition.
- NB :** Ignoring this protocol can result in the failure of the surgery and in the loss of the medical device.  
**NB2 :** We strongly recommend **not to use** automatic washing machines.

Implants Diffusion International - 23/25 rue Emile Zola 93100 Montreuil France  
+33 (0)1.48.70.70.48 [www.idi-dental.fr](http://www.idi-dental.fr)