

Brand	Reference	Catalog Length (in mm)	Absolute Total Length (in mm)	Collar Height (in mm)	Apical Diameter (in mm)	Absolute Maximum Diameter (in mm)	3D Model File
ID Max	IM+0640	6.00	7.20	0.00	4.10	4.40	ID\IM+0640.impl
ID Max	IM+0650	6.00	7.20	0.00	4.60	4.90	ID\IM+0650.impl
ID Max	IM+0840	8.00	8.00	0.00	4.10	4.40	ID\IM+0840.impl
ID Max	IM+0850	8.00	8.00	0.00	4.60	4.90	ID\IM+0850.impl
ID Max	IM+0860	8.00	8.00	0.00	6.00	6.60	ID\IM+0860.impl
ID Max	IM+1040	10.00	10.00	0.00	4.10	4.40	ID\IM+1040.impl
ID Max	IM+1050	10.00	10.00	0.00	4.60	4.90	ID\IM+1050.impl
ID Max	IM+1060	10.00	10.00	0.00	6.00	6.60	ID\IM+1060.impl
ID Max	IM+1240	12.00	12.00	0.00	4.10	4.40	ID\IM+1240.impl
ID Max	IM+1250	12.00	12.00	0.00	4.60	4.90	ID\IM+1250.impl
ID Max	IM+1260	12.00	12.00	0.00	6.00	6.60	ID\IM+1260.impl
ID Max	IM+1540	15.00	15.00	0.00	4.10	4.40	ID\IM+1540.impl
ID Max	IM+1550	15.00	15.00	0.00	4.60	4.90	ID\IM+1550.impl
ID Max	IM+1560	15.00	15.00	0.00	6.00	6.60	ID\IM+1560.impl
IDAII	IDA1042	10.00	10.00	0.00	4.20	4.20	ID\IDA1042.impl
IDAII	IDA1052	10.00	10.00	0.00	5.20	5.20	ID\IDA1052.impl
IDAII	IDA1242	12.00	12.00	0.00	4.20	4.20	ID\IDA1242.impl
IDAII	IDA1252	12.00	12.00	0.00	5.20	5.20	ID\IDA1252.impl
IDAII	IDA1542	15.00	15.00	0.00	4.20	4.20	ID\IDA1542.impl
IDAII	IDA1552	15.00	15.00	0.00	5.20	5.20	ID\IDA1552.impl
IDAII	IDA1842	18.00	18.00	0.00	4.20	4.20	ID\IDA1842.impl
IDAII	IDA1852	18.00	18.00	0.00	5.20	5.20	ID\IDA1852.impl
IDBio	IDB0835	8.00	8.3	0.00	3.50	3.70	ID\IDB0835.impl
IDBio	IDB0840	8.00	8.3	0.00	4.00	4.20	ID\IDB0840.impl
IDBio	IDB0850	8.00	8.3	0.00	4.90	5.20	ID\IDB0850.impl
IDBio	IDB1035	10.00	10.00	0.00	3.50	3.70	ID\IDB1035.impl
IDBio	IDB1040	10.00	10.00	0.00	4.00	4.20	ID\IDB1040.impl
IDBio	IDB1050	10.00	10.00	0.00	4.90	5.20	ID\IDB1050.impl
IDBio	IDB1235	12.00	12.00	0.00	3.50	3.70	ID\IDB1235.impl
IDBio	IDB1240	12.00	12.00	0.00	4.00	4.20	ID\IDB1240.impl
IDBio	IDB1250	12.00	12.00	0.00	4.90	5.20	ID\IDB1250.impl
IDBio	IDB1535	15.00	15.00	0.00	3.50	3.70	ID\IDB1535.impl
IDBio	IDB1540	15.00	15.00	0.00	4.00	4.20	ID\IDB1540.impl
IDBio	IDB1550	15.00	15.00	0.00	4.90	5.20	ID\IDB1550.impl
IDCam M	IDCM0842	8.00	8.50	0.00	4.20	4.20	ID\IDCM0842.impl
IDCam M	IDCM0852	8.00	8.50	0.00	5.20	5.20	ID\IDCM0852.impl
IDCam M	IDCM1042	10.00	9.40	0.00	4.20	4.20	ID\IDCM1042.impl
IDCam M	IDCM1052	10.00	9.40	0.00	5.20	5.20	ID\IDCM1052.impl
IDCam M	IDCM1242	12.00	11.40	0.00	4.20	4.20	ID\IDCM1242.impl
IDCam M	IDCM1252	12.00	11.40	0.00	5.20	5.20	ID\IDCM1252.impl
IDCam M	IDCM1542	15.00	14.40	0.00	4.20	4.20	ID\IDCM1542.impl
IDCam M	IDCM1552	15.00	14.40	0.00	5.20	5.20	ID\IDCM1552.impl
IDCam S	IDCS0842	8.00	8.50	0.00	4.20	4.20	ID\IDCS0842.impl
IDCam S	IDCS0852	8.00	8.50	0.00	5.20	5.20	ID\IDCS0852.impl
IDCam S	IDCS1035	10.00	9.40	0.00	3.50	3.70	ID\IDCS1035.impl
IDCam S	IDCS1042	10.00	9.40	0.00	4.20	4.20	ID\IDCS1042.impl
IDCam S	IDCS1052	10.00	9.40	0.00	5.20	5.20	ID\IDCS1052.impl
IDCam S	IDCS1235	12.00	11.40	0.00	3.50	3.70	ID\IDCS1235.impl
IDCam S	IDCS1242	12.00	11.40	0.00	4.20	4.20	ID\IDCS1242.impl
IDCam S	IDCS1252	12.00	11.40	0.00	5.20	5.20	ID\IDCS1252.impl
IDCam S	IDCS1535	15.00	14.40	0.00	3.50	3.70	ID\IDCS1535.impl
IDCam S	IDCS1542	15.00	14.40	0.00	4.20	4.20	ID\IDCS1542.impl
IDCam S	IDCS1552	15.00	14.40	0.00	5.20	5.20	ID\IDCS1552.impl
IDOne T	IDOT0840	8.00	8.00	0.00	4.10	4.40	ID\IDOT0840.impl
IDOne T	IDOT0850	8.00	8.00	0.00	4.60	4.90	ID\IDOT0850.impl

IDOne T	IDOT1040	10.00	10.00	0.00	4.10	4.40	ID\IDOT1040.impl
IDOne T	IDOT1050	10.00	10.00	0.00	4.60	4.90	ID\IDOT1050.impl
IDOne T	IDOT1240	12.00	12.00	0.00	4.10	4.40	ID\IDOT1240.impl
IDOne T	IDOT1250	12.00	12.00	0.00	4.60	4.90	ID\IDOT1250.impl
IDOne T	IDOT1540	15.00	15.00	0.00	4.10	4.40	ID\IDOT1540.impl
IDOne T	IDOT1550	15.00	15.00	0.00	4.60	4.90	ID\IDOT1550.impl