



IMPLANTS DIFFUSION INTERNATIONAL
30 years OF IMPLANTOLOGY

TURBODrill®




« INTENSE » IRRIGATION






Drill differently !

The **TURBOdrill**[®] drill is an innovation of the **IDI** company – **IMPLANTS** DIFFUSION INTERNATIONAL. It offers a new drilling approach thanks to its secure bladed stopper.




Control each drilling parameters :
safety, irrigation,
temperature & speed.

- Safety** - Titanium stopper
- Innovation** - Bladed stopper
- Efficiency** - 4 serrated blades



An abundant irrigation down into the socket

The «hollow» stopper concentrates and guides the irrigation liquid along the blades of the drill. **The more the rotation speed is high, the more the liquid is driven down into the socket.**



Optimized irrigation.

Drilling until 1200 rpm **without bone heating**. The liquid is driven by the blades down into the socket for an **optimal cooling**.

Cylindro-tapered drill,
available in 4 lengths :
8, 10, 12 et 15 mm
and 3 diameters :
3,6, 4,2 et 5,2 mm.



IMPLANTS DIFFUSION INTERNATIONAL

French designer and manufacturer of dental implants

Since 1987, **IDI** is offering innovative and quality products, **all developed and manufactured in France**, and meeting the national and international standards and certifications.

IDI – **IMPLANTS** DIFFUSION INTERNATIONAL
23/25 rue Émile Zola - 93100 Montreuil - France
Tél. : +33 (0)1 48 70 70 48 – Fax : +33 (0)1 48 70 44 58
Mail : export@idisystem.fr

www.idi-dental.com